

# Ventilated capsule-type flow compensating digital perspirometer

# SMN-1000SR

## Portable and easy to carry.

- The A5-sized case makes the SMN-1000SR easy to carry and use.
- It can be connected to a computer via a single cable.
- It supports data storage on memory cards, and the built-in battery enables usage without connecting to an external power source or device.
- It supports wireless communication using a Bluetooth LE-compatible PC or smartphone. (SMN-1000SR only)



SKINOS  
Technology  
& Products

## Uses the flow compensation method.\*

- It can stably measure heavy sweating.
  - It is capable of high-responsivity measurements.
- \*Patented technology


## What is a ventilated capsule-type perspirometer?

A ventilated capsule-type perspirometer is a simple perspiration meter that does not require dry air as it uses the air in the measurement environment as a reference. Air is supplied to the capsule that covers the skin, and the amount of sweat is measured from the humidity difference between the air before and after it passes through the skin.


## What is the flow compensation method?

When there is a small amount of sweating, the air flow rate is reduced to increase sensitivity, and when there is a large amount of sweating, the airflow rate is increased to increase the dynamic range.

Changing the airflow rate in this manner, it is possible to stably measure from small to large amounts of sweating. The output value is adjusted to be linear with respect to the amount of sweating (the amount of water evaporated from the skin).



Use example: measurement of chest sweat rate and behavior in summer



Ventilated capsule-type flow compensating digital perspirometer


# SMN-1000SR

Science of sweat  
**SKINOS**



## Specifications

Product name	SMN-1000S	SMN-1000SR
<b>Major specifications</b>		
Operating mechanism	Ventilated capsule-type differential method	
Number of channels	1Ch	
Flow rate	Flow rate compensation method (300-1,000mL)	
Analog output	Voltage output $\pm 5 V$ *1V = 1 mg/min. The voltage [V] recorded by the recorder can be converted directly into the amount of sweat [mg/(cm <sup>2</sup> · min)].	
Digital output	USB	
Wireless output	—	Bluetooth LE
Display	Measurement value display	
Memory function	Memory card	
Power supply	6 V (for battery charging)	
Built-in battery	2000 mAh (continuous operation time: 4 hours or so)	
Accessories	<ul style="list-style-type: none"> <li>• Main unit</li> <li>• Data transmission cable</li> <li>• Dedicated memory card</li> <li>• Recording analysis software</li> <li>• AC adapter charger</li> </ul>	
Recommended use conditions	Temperature: 23.5°C Humidity: 60% RH	
<b>Properties</b>		
Start-up time	~1 s	
Measurement range	0-3 mg/(cm <sup>2</sup> /min)	
Accuracy (error rate)	$\pm 10\%$ (under recommended conditions, range 0-1 mg/(cm <sup>2</sup> /min))	



\*This device is designed for research applications.

Science of sweat  
**SKINOS**

**SKINOS Co., Ltd.**

2-16-24 Fumiiri, Ueda City, Nagano Prefecture, 386-0017 Japan  
Shinshu University Open Venture Innovation Center, Room 107  
TEL : 0268-75-9071 FAX : 0268-75-9072 info@skinos.co.jp https://www.skinos.co.jp

**SHINSHU UNIVERSITY**

Startup company originating from Shinshu University  
JAPAN

